

Generic External Condom Catheter Parts List

This list provides the links to the individual components to create your own “Stadium Pal” type external condom catheter setup. The IPA receives no remuneration from any product vendor.

The urine collector leg bag

<https://www.vitalitymedical.com/conveen-security-leg-bag-with-clamp-or-lever-tap-outlet.html>

High-quality leg straps to secure the leg bag

<https://www.vitalitymedical.com/coloplast-conveen-security-fabric-leg-bag-straps.html>

The condom with drainage tip

<https://www.vitalitymedical.com/coloplast-freedom-clear-catheter.html>

Skin Prep – use in advance of applying the condom to allow for easier removal after use

<https://www.vitalitymedical.com/smith-and-nephew-skin-prep-protective-dressing.html>

Uni-Solve adhesive remover – for after use

<https://www.vitalitymedical.com/smith-and-nephew-unisolve-adhesive-remover.html>

Syringe to flush the bag easily--inject mix of a few ounces of water with a tablespoon of bleach

<https://www.vitalitymedical.com/amsino-60-ml-thumb-control-ring-syringe-as016-as015.html>

Compiled by Dan Rocker, September 2025